


Custom Walker Adjustment

With your arm extended down, handle should be at wrist height.
With your hand placed on the hand-grip, your elbow should be flexed 20° - 30°.



Beware of:

- Electrical cords in way
- Liquids on floor
- Throw rugs
- Toys or debris
- Furniture
- Doors and door stops
- Gaps in floors

Other Precautions:

- 1 Use wheels on front legs only or on all four legs with appropriate rear wheelbrakes.
- 2 Do not push one side of the walker ahead and then other side without lifting the walker first.
- 3 Use Merits accessories only. Use of other manufacturers accessories can result in improper fit and risk fall and/or injury.
- 4 Do not use on stairs or escalators.

Lifetime Limited Warranty

The product is warranted to be free of defects in materials and workmanship for the lifetime of the original consumer purchaser.

This device was built to exacting standards and carefully inspected prior to shipment. This Lifetime Limited Warranty is an expression of our confidence in the materials and workmanship of our products and our assurance to the consumer of years of dependable service.

In the event of a defect covered by this warranty, we will, at our option, repair or replace the device.

This warranty does not cover device failure due to owner misuse or negligence, or normal wear and tear. The warranty does not extend to non-durable components, such as rubber accessories, casters, and grips, which are subject to normal wear and need periodic replacement.

If you have a question about your product or this warranty, please contact an authorized dealer.



MERITS HEALTH PRODUCTS, INC.

P.O. BOX 150356
Cape Coral, FL 33915
1-800-963-7487
1-800-(9-MERITS)

www.meritshealth.com



Folding Mechanism



Folding Walker

Instructions & Warranty Information

For use with all Two Button Walkers.

Instructions for easy patient transfer



merits

Wheelchair to Walker Transfer



*M*aking certain wheels on the wheelchair are locked, lift your body by pushing down on the armrests and carefully transfer one hand at a time to the walker.




Wheelchair to Walker Transfer

(with lower extremity injury to right leg or foot)



*M*aking certain wheels on the wheelchair are locked, lift your body by pushing down on the armrests and carefully transferring first the left hand, then the right hand to the walker.
(If injury to left leg or foot transfer the right hand first, followed by left.)



Three-Point Gait (for walking with walker)

*A*lign the middle of your foot with the back legs of the walker...



*L*ift the walker up and place it comfortably in front of you, making sure all four legs end up on the floor...



*T*hen carefully walk forward, one step at a time.

