

LIGHTWEIGHT PORTABLE SCOOTER



Seat Color  
Black Vinyl

Body Colors  
Blue  
Red

**New!**

# Ride

**BUZZ** around **XL**

3 wheel  
GB-116

Golden proudly introduces the all-new Buzzaround XL! Featuring a 300 lb. weight capacity, a generous seat size and completely wireless disassembly, the Buzzaround XL is the latest generation of the popular Buzzaround series. Available in Red and Blue, the XL comes standard with a Delta tiller and a battery charging port is conveniently located on the tiller under the control panel.



Rise. Roll. Rest. Ride.

**GOLDEN**  
technologies

YOUR LIFE IN MOTION



### 3 wheel Specifications

<b>Model</b>	<b>GB116</b>
Medicare Code	K0800 (pending)
Weight Capacity	300 lbs
Drive Wheels	Rear
Maximum Speed	4.5 mph
Operating Range <sup>1</sup>	5.5 miles
Ground Clearance to center deck	2.25"
Ground Clearance to transaxle/motor	1.5"
Turning Radius	37"
Type Batteries	2-12AH
Colors	Red, Blue
Freewheel Mode	Yes
Electronic Speed Control	Yes
Electro-Mechanical Brakes	Yes
Charger	Off Board
Controller	50 Amp Dynamic R Series
<b>Dimensions:</b>	
Length	40.5"
Width	21.5"
Height (ground to top of back of seat)	34.75" - 36.75"
Ground to Top of Seat	22" - 24"
Ground to Top of Deck	4.125"
Deck to Top of Seat	17.25" - 19.25"
Number of Seat Height Positions	4
Size of Increments	0.5"
Front Axle to Rear Axle	30"
<b>Weight of Unit: (assembled)</b>	
Front Half	33 lbs
Rear Half	25 lbs
Seat w/Arms	21 lbs
Batteries	21 lbs/pack
<b>Tires:</b>	
Front	8"
Rear	8"
Rear Anti-Tippers	2"
<b>Standard Seat:</b>	
Slide Seat	Stadium Style
Back Height - no headrest	N/A
Back Height with headrest	16"
Width x Depth	17" X 16"
Color	Black Vinyl
<b>Optional Seats:</b>	
Back Height - no headrest	Stadium Style
Back Height with headrest	16"
Width x Depth	18" x 16", 20" x 18"
Color	Black Vinyl
<b>Warranty:</b>	
Frame	Lifetime
Drive Train	2 Years
Electronics <sup>2</sup>	2 Years

### Wireless disassembly in 3 easy steps!



**Step 1:** Pop off the seat!



**Step 2:** Pull out the battery pack!



**Step 3:** Pull up on the drivetrain release lever to separate the front and rear frame sections!



Disassembles into 5 pieces!



A battery charging port is conveniently located on the tiller under the control panel.

### Features

- Exclusive interlocking teeth for fully connectorless assembly/disassembly
- 300 lb. weight capacity
- Auto grade cut pile carpet
- Tiller mounted basket

<sup>1</sup> Battery range will vary due to rider weight, drive surface, terrain and battery size.  
<sup>2</sup> Electronics warranty excludes batteries.

Literature is current at time of printing. Golden Technologies reserves the right to make changes to the product or literature at any time.

