


VALUE SERIES (cost-effective)

This series combines ultra-soft lumbar pillow support with our unique "grip-rite" arm design for easy transfer.



*Std. on PR200, PR355M & PR355L



Monarch Plus Md.
PR-359M
shown in Cornflower



Monarch Plus Lg.
PR-359L
shown in Rosewood



Monarch Md.
PR-355M
shown in Autumn




Monarch Lg.
PR-355L
shown in Spring



Capri
PR-200
shown in Chestnut

TRADITIONAL SERIES (universal designs)



*Std. on PR401M, PR401L & PR747

The Traditional series is all about looks. Choose from colonial or overstuffed biscuit back classics.



Cambridge
PR-401
shown in Sterling




Cambridge
PR-401L
shown in Copper



Williamsburg
PR-747
shown in Island

SIGNATURE SERIES (multi-functional)

Want a little something extra from your lift chair? This selection has it all, and... even a little more.



Winston
PR-410
shown in Saddle

The overstuffed, biscuit back and plush, rounded armrests offer generous proportions for roomy relaxation.



Regal
PR-751

Ease into this luxurious 3-pillow lift chair with convenient arm storage and foldable tray.



Space Saver
PR-906

This unique frame design allows for chair placement to be only 6" from the wall, the perfect chair for compact living spaces!





Royal
PR-752

The 3-pillow waterfall back and full luxury chaise pad is a classic style offering full support.



Transfer Chair
PR-451

Unique flip-up arm design allows for easy transfer into and out of a power chair.



PRODUCT SPECIFICATIONS

Model	Weight Capacity	Overall Width	Overall Height	Seat to Top of Back	Width Between Arms	Seat Depth	Floor to Top of Seat	Extension, While Reclined from Back of Seat to End of Footrest	Floor to Top of Footrest	Heat & Massage Option Available	Footrest Extension Option Available	Straight Lift Option Available*
Traditional												
Cambridge	PR-401	375 lbs.	33"	41.5"	28"	20"	20"	19"	41.5"	22"	Yes	Yes
Cambridge Lg.	PR-401L	375 lbs.	35"	43.5"	30"	22"	22"	21.5"	44"	23"	Yes	Yes
Williamsburg	PR-747	375 lbs.	32.5"	44"	28"	21"	21.5"	19.5"	42"	23"	Yes	Yes
Signature												
Winston	R-410	375 lbs.	43"	46"	29"	21.5"	21"	21"	45"	23.5"	Yes	No
Transfer	PR-451	375 lbs.	33"	42.5"	29"	20"	20"	19"	40.5"	23"	Yes	Std.
Regal	PR-751TY	375 lbs.	36"	44"	28"	22"	21.5"	20"	42"	24"	Yes	No
Royal	PR-752	375 lbs.	32.5"	44"	28"	21"	21"	19.5"	42"	23.5"	Yes	Yes
Space Saver	PR-906	250 lbs.	35.5"	43"	27"	18"	22"	22"	44"	21.5"	Yes	No
MaxiComfort												
Relaxer Md.	PR-756MC	375 lbs.	36.5"	44"	29"	20"	21"	19"	43"	23"	Yes	Yes
Relaxer Lg.	PR-756L	375 lbs.	39"	46"	30"	22"	22"	20.5"	44"	24"	Yes	Yes
Windsor	PR-516	375 lbs.	39"	45"	28"	21"	20"	21"	44"	24"	Yes	No
Cloud	PR-510	375 lbs.	38"	45"	29"	20"	21.5"	20"	44"	24.5"	Yes	Yes
Comforter JP	PR-505JP	300 lbs.	30.5	39"	22.5"	19"	18.5"	18"	38"	24.5"	Yes	No
Comforter Sm.	PR-505S	300 lbs.	31.5	40"	25"	20"	19.5"	19"	39"	23"	Yes	Yes
Comforter Md.	PR-505M	375 lbs.	32.5	42"	27"	21"	21"	19.5"	42"	23"	Yes	Yes
Comforter Lg.	PR-505L	375 lbs.	34.5	44"	28.5"	23"	22"	20.5"	43"	23"	Yes	Yes
Comforter												
Junior Petite	R-501JP	300 lbs.	30.5"	39"	22.5"	19"	18.5"	18"	38"	22"	Yes	No
Small	PR-501S	300 lbs.	31.5"	40"	25"	20"	19.5"	19"	39"	23"	Yes	Yes
Medium	PR-501M	375 lbs.	32.5"	42"	27"	21"	21"	19.5"	42"	23"	Yes	Yes
Large	PR-501L	375 lbs.	34.5"	44"	28.5"	23"	22"	20.5"	43"	23"	Yes	Yes
Tall	PR-501T	375 lbs.	34.5"	47"	30"	23"	24"	22.5"	45"	23"	Yes	Yes
Comforter Wides												
Small X-wide	PR-501S-23	375 lbs.	34.5"	40"	25"	23"	19.5"	18"	39"	23"	Yes	Yes*
Medium-26 Dbl.	PR-501M-26D	500 lbs.	38.5"	42"	27"	27.5"	21.5"	19.5"	42"	23"	Yes	Yes*
Large-26 Dbl.	PR-501L-26D	500 lbs.	38.5"	44"	28.5"	27.5"	22.5"	20.5"	43"	23"	Yes	Yes*
Tall-28 Dbl.	PR-501T-28D	500 lbs.	41"	48"	30"	29.5"	24"	22.5"	45"	23"	Yes	Yes*
Super 33	PR-502	700 lbs.	45"	47"	29"	33"	22"	21"	45"	21"	Yes	No
Value												
Capri	PR-200	325 lbs.	31"	41"	27"	20"	19"	18"	40"	23"	No	No
Monarch Md.	PR-355M	375 lbs.	32"	42"	28"	21"	21.5"	18.5"	43"	23"	Yes	No
Monarch Lg.	PR-355L	375 lbs.	34"	44"	30"	22"	22"	18.5"	43"	22"	Yes	No
Monarch Plus Md.	PR-359M	375 lbs.	32"	41.5"	27.5"	20"	21"	19.5"	43"	23"	Yes	Yes
Monarch Plus Lg.	PR-359L	375 lbs.	34"	43"	29.5"	23"	22"	20.5"	43"	23"	Yes	Yes

*Must be set at factory.



Visit Your Authorized Golden Dealer:



BR-LIFT-121812



Golden Lift Chairs

AMERICA'S #1 SELLING BRAND



Our American made lift chairs come in every shape and size for the perfect fit.

Our premium fabrics match nearly every décor, and our patented technologies give you the most comfortable lift and recline experience.



www.goldentech.com

MAXICOMFORT® SERIES

(maximum comfort)




These are the most advanced lift and recline chairs available anywhere in the world!



MaxiComfort® chairs transport your body into a blissful, relaxed state. Choose from optimal comfort positions that promote health, wellness, and relaxation for both your body and your mind.

MAXICOMFORTER FABRICS



MaxiComforter Jr. Petite
PR-505JP
shown in Cabernet

MaxiComforter Small
PR-505S
shown in Palomino


MaxiComforter Medium
PR-505M
shown in Admiral

MaxiComforter Large
PR-505L
shown in Evergreen

COMFORTER SERIES

(a perfect fit)

The Comforter series fits people who are 5'0" to 6'8" tall with weight capacities from 300 to 375 lbs. These chairs give the perfect sizes for the perfect fit.



Jr. Petite
PR-501JP
shown in Admiral
For heights 5'0" and under

Small
PR-501S
shown in Pearl
For heights 5'0" to 5'3"

Medium
PR-501M
shown in Palomino
For heights 5'4" to 5'10"

Large
PR-501L
shown in Evergreen
For heights 5'11" to 6'2"


Tall
PR-501T
shown in Cabernet
For heights 6'3" to 6'8"



MaxiComfort Cloud
PR-510
shown in Sterling

When you sit in the Cloud, you feel cushion and support from every angle. Get surrounded with first class comfort.


CLOUD FABRICS



Windsor
PR-516
shown in Coffee Bean Brisa

The Windsor offers contemporary styling with a plush, two pillow back and overstuffed flared armrests.

WINDSOR FABRICS



Relaxer
PR-756MC
shown in Vino


RELAXER FABRICS



Relaxer
PR-756L
shown in Saddle

COMFORTER WIDES

The Comforter Wides are perfect for those who need more space between the armrests. Constructed using the highest quality components for reliable, steady lifting and reclining, these chairs offer weight capacities from 375 to 700 lbs.



Small 23
PR-501S-23
shown in Admiral
For heights 5'0" to 5'3"

Medium 26
PR-501M-26
shown in Palomino
For heights 5'4" to 5'10"
(Dual Motor)

Large 26
PR-501L-26
shown in Cabernet
For heights 5'11" to 6'2"
(Dual Motor)

Tall 28
PR-501T-28 shown in Evergreen
For heights 6'3" to 6'8"
(Dual Motor)

Super 33
PR-502
shown in Admiral
(Triple Motor)



Travel effortlessly and smoothly through your favorite positions. With pre-programmed and programmable controls, you can lift, recline, watch TV, sleep, or enjoy the wellness benefits of the zero-gravity and Trendelenburg positions, all at the touch of the AutoDrive buttons.