

R E L A X E R S E R I E S


Rise




Relaxer PR-756L



WINDSOR PR-516 FABRICS



RELAXER PR-756 FABRICS



Windsor PR-516



Relaxer PR-756MC

Relaxer Series

MaxiComfort is now available on more models than ever before! Golden's patented MaxiComfort Ultimate Recline Technology is standard on the Relaxer Medium and Relaxer Large, as well as the Windsor lift chair. Every MaxiComfort chair comes standard with our patented MaxiComfort Ultimate Recline Technology for correct spinal alignment, back pressure relief and improved circulation.



Rise. Roll. Rest. Ride.



YOUR LIFE IN MOTION

RELAXER SERIES



	Relaxer	Relaxer	Windsor
Model Number	PR-756MC	PR-756L	PR-516
Size	Medium	Large	Medium
Weight Capacity	375 LB	375 LB	375 LB
Recline Positions	MaxiComfort	MaxiComfort	MaxiComfort
Standard Fabrics	4	4	4
Fabric Options	32	32	Limited
Chaise Seat	YES	YES	YES
Back Type	4 Pillow	4 Pillow	2 Pillow
Head Rest Cover	YES	YES	YES
Arm Covers	NO	NO	NO
Battery Back Up	YES	YES	YES
Overall Width	36.5"	39"	39"
Overall Height	44"	46"	45"
Width Between Arms	20"	22"	21"
Floor to Top of Seat	19"	20.5"	21"
Seat Depth	21"	22"	20"
Seat to Top of Back	29"	30"	28"
Seat to Top of Arm	8"	8"	5"
Floor to Back Edge of Seat Raised	26.5"	28"	28"
Floor to Front Edge of Seat Raised	22"	23"	24.5"
Seat Angle Raised	15°	15°	16°
Floor to Top of Foot Rest Reclined	23"	24"	24"
Floor to Top of Back Reclined	21"	22"	25"
Back Angle Closed	100°	100°	95°
Back Angle Reclined	174°	174°	165°
Foot Rest Angle Closed	90°	90°	90°
Foot Rest Angle Opened	180°	180°	180°
Distance Required from Wall Reclined	25"	27"	28"
Extension Reclined from Back of Seat to End of Foot Rest	43"	44"	44"
Life Time Warranty, Chair Frame, Lift Frame, Recline Mechanism	YES	YES	YES
Warranty on Electrical Parts	3 Yrs.	3 Yrs.	3 Yrs.
Pro Rated Years 4 Thru 7	YES	YES	YES
Warranty on Mechanical Labor	3 Yrs.	3 Yrs.	3 Yrs.
Warranty on Electrical Labor	1 Yr.	1 Yr.	1 Yr.
Removable Back Option Available	NO	NO	NO
Smart Tek System Standard	YES	YES	YES
Customizable	NO	NO	NO
Left Hand Control Available	YES	YES	YES
Extra Pocket Option Available	YES	YES	YES
Heat and Massage Option Available	YES	YES	YES
Foot Rest Extension Option Available	YES	YES	NO
Straight Lift Option Available	NO	NO	NO
User Height (suggested)	5'4"-5'10"	5'11"-6'2"	5'4"-5'10"

Sizes may vary depending on fabric, filling material, upholstery and carpet thickness. All measurements taken on concrete floor using levels and metal rulers.



501 M shown

Golden Craftsmanship

Each luxury lift & recline chair is hand crafted in our facility using decades of motion furniture experience.



501 M shown

Added Comfort and Safety

Designed to perform year after year using advanced technology backed by the best warranty.



501 M shown

Smart Engineering

Smooth, quiet lifting system offers a secure lift for the best performance.



Smart Tek

Our luxurious lift and recline chairs are equipped with electronics to help assure a secure lift every time you need it.



Literature is current at time of printing. Golden Technologies reserves the right to make changes to the product or literature at any time.



Complies with UL 73 and CSA C22.2 No. 68



Listing # E112417



GOLDEN technologies