

Backpacker[®] MV

Base weight lift 191 lbs.¹

Scooter/power chair overall weight 350 lbs.

Maximum product width 27"

Maximum travel distance from the ground 37"

Motor Sealed 12 volt DC

Overall platform length 44" (41" usable surface)

Maximum product wheelbase 43"

Warranty 3-year transferable

Optional battery pack 12 volt, on-board

Optional charger 110 volt

¹ Weight does not include any accessories

² Vehicle dependent

Can easily lift the following current Pride[®] Power Chairs²

- Jazzy[®] 600
- Jazzy[®] 614/614 HD
- Jazzy[®] 1113 ATS
- Jazzy Select[®]
- Jazzy Select[®] 6
- Jazzy Select[®] 6 Ultra
- Jazzy Select[®] Traveller
- Jazzy Select[®] Elite
- Jet[™] 3 Ultra
- Go-Chair
- Z-Chair

Can easily lift the following current Quantum[®] Power Chairs²

- Quantum[®] 600/600 XL/640
- Quantum[®] 610
- Quantum[®] 614
- Quantum[®] 1103 Ultra
- Quantum[®] 1121
- Quantum[®] 6400
- Quantum[®] 6000Z/6400Z/6000Z HD
- Quantum[®] R4000/R4400

Can easily lift the following current Pride[®] Scooters²

- Celebrity[®] X 3 & 4
- Go-Go[®] Ultra X 3 & 4
- Go-Go Elite Traveller[®]/LR 3 & 4
- Go-Go Elite Traveller[®] Plus 3 & 4
- Go-Go Elite Traveller[®] Plus HD 3 & 4
- Legend[®] 3 & 4
- Maxima 3 & 4
- Revo[®] 3 & 4
- Victory[®] 9 3 & 4
- Victory[®] 10 3 & 4
- Victory[®] 9 PS

¹ Weight does not include any accessories

² Vehicle dependent

Battery Pack

NEW, OPTIONAL MULTI-PURPOSE BATTERY BOX



- Industry first
- Industry exclusive
- Charges the mobility device from the lift as you drive

Backpacker[®] MV



Backpacker[®] MV in lowest position.

Silver Star[®] Lifts & Ramps

www.silverstarmobility.com • www.pridemobility.com
(800) 800-8586 (U.S.) • (888) 570-1113 (Canada)

Backpacker[®] MV

Now you can take your mobility on the road easily in your Sport Utility Vehicle with the Backpacker[®] MV from Silver Star Lifts and Ramps. It's the ultimate solution to lifting your powered mobility device into your vehicle, engineered specifically to be a perfect fit for a wide variety of popular vehicles including SUVs. Like all Silver Star Lifts, the Backpacker MV features all-steel construction and a protective finish, along with an industrial rated 12 volt motor and battery pack, and easy-to-use controls to combine rugged durability with safe, simple and reliable operation.

Note: The Backpacker[®] MV can lift all of the listed power chairs with contoured, Synergy[®] or Specialty seating. Please call Customer Service when requesting information about lifting with a power positioning system. Can fit into most SUV's and mini-vans with little or no vehicle modification, full-size vans require modification.

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.



Backpacker[®] MV fully loaded.

- One piece, fully assembled drive-on/drive-off platform
- Removable, front drive panel for easier serviceability
- New, adjustable wheel chocks on the platform surface
- Single piece horizontal drive rail: Minimum welding points for noise reduction, increased rigidity, and more fluid movement
- New, angled, low-profile platform
- Easy-to-use 2 button hand control
- Can lift a wide range of Pride[®] Scooters and Power Chairs and other manufacturers' mobility devices

- 350 lbs. lifting capacity
- 3-year transferable warranty
- Fits into SUVs with a 31" opening height
- Manual override
- Improved micro switch protection
- False floor available for specified vehicles through special order to preserve vehicle aesthetics
- Ratcheting straps now available as an option