

**PRIDE
TOUGH**

Portable Ramps



Single-Fold Ramps



Solid Ramps

Portable Ramps

Ramps

Portable Ramps Chart

° = degree of slope
% = percent grade

Note: 6', 12' Solid Ramps and all 34" wide ramps are not available with Grip tape

Ramp Length	6' (15 cm)	12' (30 cm)	1.5' (46 cm)	2' (61 cm)	3' (91 cm)	4' (122 cm)	5' (152 cm)	6' (183 cm)	Rise
Solid Ramps	19.5° (35%)	12.8° (23%)	16.1° (29%)	14.5° (26%)	17° (30%)	14.5° (26%)	12.8° (23%)	16.1° (29%)	20" (51 cm)
	14.5° (26%)	9.6° (17%)	12.8° (23%)	9.6° (17%)	14.5° (26%)	11.2° (20%)	12.8° (23%)	14.5° (26%)	18" (46 cm)
	9.6° (17%)	6.4° (11%)	9.6° (17%)	6.4° (11%)	9.6° (17%)	8.4° (15%)	9.6° (17%)	16" (41 cm)	16" (41 cm)
	4.8° (8%)	3.2° (6%)	6.4° (11%)	4.8° (8%)	6.4° (11%)	7.2° (13%)	7.2° (13%)	14" (36 cm)	14" (36 cm)
	-	-	4.8° (8%)	3.2° (6%)	4.8° (8%)	6.4° (11%)	6.4° (11%)	12" (30 cm)	12" (30 cm)
	-	-	3.2° (6%)	2.4° (4%)	3.2° (6%)	4.8° (8%)	4.8° (8%)	10" (25 cm)	10" (25 cm)
	-	-	2.4° (4%)	1.9° (3%)	2.4° (4%)	3.2° (6%)	3.2° (6%)	8" (20 cm)	8" (20 cm)
	-	-	1.9° (3%)	1.4° (2%)	1.9° (3%)	2.4° (4%)	2.4° (4%)	7" (18 cm)	7" (18 cm)
	-	-	1.4° (2%)	1.0° (2%)	1.4° (2%)	1.9° (3%)	1.9° (3%)	6" (15 cm)	6" (15 cm)
	-	-	1.0° (2%)	0.7° (1%)	1.0° (2%)	1.4° (2%)	1.4° (2%)	5" (13 cm)	5" (13 cm)
Single-Fold Ramps	16.1° (29%)	12.8° (23%)	17° (30%)	14.5° (26%)	16.1° (29%)	14.5° (26%)	12.8° (23%)	16.1° (29%)	18" (46 cm)
	14.5° (26%)	9.6° (17%)	17° (30%)	14.5° (26%)	16.1° (29%)	14.5° (26%)	12.8° (23%)	14.5° (26%)	16" (41 cm)
	12.8° (23%)	8.4° (15%)	14.5° (26%)	11.2° (20%)	14.5° (26%)	11.2° (20%)	12.8° (23%)	14.5° (26%)	14" (36 cm)
	10.8° (19%)	8.4° (15%)	12.8° (23%)	9.6° (17%)	12.8° (23%)	10.8° (19%)	12.8° (23%)	12" (30 cm)	12" (30 cm)
	9.6° (17%)	7.2° (13%)	10.8° (19%)	8.4° (15%)	10.8° (19%)	9.6° (17%)	10.8° (19%)	10" (25 cm)	10" (25 cm)
	8.4° (15%)	6.4° (11%)	9.6° (17%)	7.2° (13%)	9.6° (17%)	8.4° (15%)	9.6° (17%)	8" (20 cm)	8" (20 cm)
	7.2° (13%)	5.7° (10%)	8.4° (15%)	6.4° (11%)	8.4° (15%)	7.2° (13%)	8.4° (15%)	7" (18 cm)	7" (18 cm)
	6.4° (11%)	4.8° (8%)	7.2° (13%)	5.7° (10%)	7.2° (13%)	6.4° (11%)	7.2° (13%)	6" (15 cm)	6" (15 cm)
	5.7° (10%)	4.8° (8%)	6.4° (11%)	4.8° (8%)	6.4° (11%)	5.7° (10%)	6.4° (11%)	5" (13 cm)	5" (13 cm)
	4.8° (8%)	3.8° (7%)	5.7° (10%)	4.8° (8%)	5.7° (10%)	4.8° (8%)	5.7° (10%)	4" (10 cm)	4" (10 cm)
3.8° (7%)	2.9° (5%)	4.8° (8%)	3.8° (7%)	4.8° (8%)	3.8° (7%)	4.8° (8%)	3" (8 cm)	3" (8 cm)	
2.9° (5%)	1.9° (3%)	3.8° (7%)	2.9° (5%)	3.8° (7%)	2.9° (5%)	3.8° (7%)	2" (5 cm)	2" (5 cm)	
1.9° (3%)	-	2.9° (5%)	1.9° (3%)	2.9° (5%)	1.9° (3%)	2.9° (5%)	1" (2.5 cm)	1" (2.5 cm)	

Pride® Lifts & Ramps

www.prideliftsandramps.com • www.pridemobility.com

St. Catharines, Ontario • (888) 570-1113

Pride®
Lifts & Ramps

PRIDE®
PERFORMANCE MOBILITY
— LIVE YOUR BEST!™ —

www.prideliftsandramps.com • www.pridemobility.com
St. Catharines, Ontario • (888) 570-1113

Single-Fold Ramps



Shown
with
grip tape

Single-Fold Ramps from Pride[®] Lifts and Ramps offer a reliable solution for transitioning curbs, thresholds and other obstacles. They are available in several different lengths and can accommodate a weight capacity of 600 lbs. Plus, lightweight Single-Fold Ramps are easily portable, making them ideal for home and travel.

Solid Ramps



Shown
with
grip tape

Solid Ramps from Pride[®] Lifts and Ramps provide long-term access for mobility devices to a variety of environments. Built with the quality you expect, Pride Solid Ramps are lightweight but strong, with a weight capacity of 600 lbs. and available in a variety of lengths.

Single-Fold Ramps

Features

- 600 lbs. weight capacity
- 3-year transferable warranty
- Lightweight & portable
- Quick delivery, nationwide distribution
- Integrated, high-traction surface
- Comprehensive product line
- Competitively priced
- 29.5" wide
- Available in: 2', 3', 4', 5' and 6'

Solid Ramps

Features

- 34" wide Solid Threshold Ramps available in 6", 12", 18" and 24"
- 600 lbs. weight capacity
- 3-year transferable warranty
- Lightweight
- Built-in traction control
- Slotted for moisture drainage
- Aluminum construction
- Side rails for safety
- Securement components for mounting purposes
- 30" wide
- Available in: 6", 12", 18", 2', 3', 4' and 5'