

**PRIDE  
TOUGH**

# Portable Ramps



Single-Fold Ramps



Solid Ramps



Multi-Fold Ramps

**Pride® Lifts & Ramps**

www.prideliftsandramps.com • www.pridemobility.com  
(800) 800-8586 (U.S.) • (888) 570-1113 (Canada)

## Ramps

# Portable Ramps

**Pride®  
Lifts & Ramps**

**PRIDE®**  
PERFORMANCE MOBILITY  
— LIVE YOUR BEST!™ —

www.prideliftsandramps.com • www.pridemobility.com  
(800) 800-8586 (U.S.) • (888) 570-1113 (Canada)

PORTBLERAMPS-BR-08-01-12

# Portable Ramps Chart

° = degree of slope  
% = percent grade

Note: 6", 12" Solid Ramps and all 34" wide ramps are not available with Grip tape

Ramp Length	Single-Fold Ramps					Multi-Fold Ramps					
	6" (15 cm)	12" (30 cm)	15" (46 cm)	2" (61 cm)	3" (91 cm)	4" (122 cm)	5" (152 cm)	6" (183 cm)	7" (213 cm)	8" (244 cm)	10" (305 cm)
19.5" (35%)	16.1" (29%)	17" (30%)	14.5" (26%)	16.1" (29%)	17" (30%)	15.5" (28%)	16.6" (30%)	17.7" (31%)	18.8" (32%)	19.9" (33%)	21" (34%)
14.5" (26%)	12.8" (23%)	14.5" (26%)	12.8" (23%)	14.5" (26%)	12.8" (23%)	14.5" (26%)	12.8" (23%)	14.5" (26%)	12.8" (23%)	14.5" (26%)	12.8" (23%)
9.6" (17%)	9.6" (17%)	9.6" (17%)	9.6" (17%)	9.6" (17%)	9.6" (17%)	9.6" (17%)	9.6" (17%)	9.6" (17%)	9.6" (17%)	9.6" (17%)	9.6" (17%)
4.8" (8%)	4.8" (8%)	4.8" (8%)	4.8" (8%)	4.8" (8%)	4.8" (8%)	4.8" (8%)	4.8" (8%)	4.8" (8%)	4.8" (8%)	4.8" (8%)	4.8" (8%)
3.2" (6%)	3.2" (6%)	3.2" (6%)	3.2" (6%)	3.2" (6%)	3.2" (6%)	3.2" (6%)	3.2" (6%)	3.2" (6%)	3.2" (6%)	3.2" (6%)	3.2" (6%)
2" (61 cm)	3" (91 cm)	4" (122 cm)	5" (152 cm)	6" (183 cm)	7" (213 cm)	8" (244 cm)	10" (305 cm)	12" (305 cm)	14" (36 cm)	16" (41 cm)	18" (46 cm)
1.9" (3%)	1.9" (3%)	1.9" (3%)	1.9" (3%)	1.9" (3%)	1.9" (3%)	1.9" (3%)	1.9" (3%)	1.9" (3%)	1.9" (3%)	1.9" (3%)	1.9" (3%)
1.5" (46 cm)	1.5" (46 cm)	1.5" (46 cm)	1.5" (46 cm)	1.5" (46 cm)	1.5" (46 cm)	1.5" (46 cm)	1.5" (46 cm)	1.5" (46 cm)	1.5" (46 cm)	1.5" (46 cm)	1.5" (46 cm)
12" (30 cm)	12" (30 cm)	12" (30 cm)	12" (30 cm)	12" (30 cm)	12" (30 cm)	12" (30 cm)	12" (30 cm)	12" (30 cm)	12" (30 cm)	12" (30 cm)	12" (30 cm)
6" (15 cm)	6" (15 cm)	6" (15 cm)	6" (15 cm)	6" (15 cm)	6" (15 cm)	6" (15 cm)	6" (15 cm)	6" (15 cm)	6" (15 cm)	6" (15 cm)	6" (15 cm)

Ramps

## Single-Fold Ramps



Shown with grip tape

Single-Fold Ramps from Pride™ Lifts and Ramps offer a reliable solution for transitioning curbs, thresholds and other obstacles. They are available in several different lengths and can accommodate a weight capacity of 600 lbs. Plus, lightweight Single-Fold Ramps are easily portable, making them ideal for home and travel.

- 600 lbs. weight capacity
- 3-year transferable warranty
- Lightweight & portable
- Quick delivery, nationwide distribution
- Integrated, high-traction surface
- Comprehensive product line
- Competitively priced
- 29.5" wide
- Available in: 2', 3', 4', 5' and 6'

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.

## Solid Ramps



Shown with grip tape

Solid Ramps from Pride™ Lifts and Ramps provide long-term access for mobility devices to a variety of environments. Built with the quality you expect, Pride Solid Ramps are lightweight but strong, with a weight capacity of 600 lbs. and available in a variety of lengths.

- 34" wide Solid Threshold Ramps available in 6", 12", 18" and 24"
- 600 lbs. weight capacity
- 3-year transferable warranty
- Lightweight
- Built-in traction control
- Slotted for moisture drainage
- Aluminum construction
- Side rails for safety
- Securement components for mounting purposes
- 30" wide
- Available in: 6", 12", 18", 2', 3', 4' and 5'

## Multi-Fold Ramps



Multi-Fold Ramps from Pride™ Lifts and Ramps make using a mobility product easy over entrances, small staircases and a variety of places in the home and outdoors. Highly durable and portable, Multi-Fold Ramps are offered in longer lengths for increased accessibility.

- 3-year transferable warranty
- Available in 6', 7', 8' at 600 lbs. and 800 lbs.\*, and 10' at 500 lbs.
- Available with or without grip tape surface
- Weight of 500 lbs. or 600 lbs. ramps with or without grip tape 6' - 32 lbs., 7' - 42 lbs., 8' - 47 lbs., 10' - 60 lbs.\*\*
- Weight of 800 lbs. ramps with grip tape 6' - 42 lbs., 7' - 47 lbs., 8' - 54 lbs.\*\*
- 29.5" useable width

\*800 lbs. ramps available only with grip tape.

\*\*Weights are approximate and may vary.