

Scooter Line Specifications

Ride



	BUZZAROUND		LITERIDER		COMPANION			AVENGER
Model	GB106	GB116	GL110	GL140	GC240	GC340	GC440	GA541
Model Name	Buzzaround Lite	Buzzaround XL	LiteRider	LiteRider	Companion	Companion	Companion	Avenger
Medicare Code	K0800	K0800	K0800	K0800	K0806	K0807	K0807	Not Classified
Weight Capacity	250 lbs	300 lbs	300 lbs	300 lbs	300 lbs	350 lbs	350 lbs	500 lbs
Drive Wheels	Rear	Rear	Rear	Rear	Rear	Rear	Rear	Rear
Maximum Speed	4 mph	4.5 mph	5 mph	5 mph	4.5 mph	4.5 mph	4.5 mph	7 mph
Operating Range ¹	8 miles	5.5 miles / 7.0 miles ⁴	9.5 miles ⁴	9.5 miles ⁴	15.5 miles	16.5 miles	15 miles	18 miles
Ground Clearance	1.75"	2.25"	3.75"	3.75"	3.75"	4.5"	4.5"	4"
Turning Radius	31"	37"	39"	48"	39"	40.5"	53.5"	67"
Type Batteries	2-12AH	2-12AH / 2-18AH	2-18AH	2-18AH	2-U1	2-U1	2-U1	2-NF22
Colors	Red, Blue, Sahara	Red, Blue	Red, Blue	Red, Blue	Red, Blue	Red, Blue	Red, Blue	Silver, Red
Freewheel Mode	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes
Electronic Speed Control	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes
Electro-Mechanical Brakes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes
Charger	Off Board	Off Board	Off Board	Off Board	Off Board	Off Board	Off Board	Off Board
Controller	45 Amp	50 Amp Dynamic R Series	90 Amp Dynamic R Series	90 Amp Dynamic R Series	90 Amp Dynamic R Series	90 Amp Dynamic R Series	90 Amp Dynamic R Series	PG 110 Amp
Dimensions:								
Length	36.5"	40.5"	41.5"	43.25"	42.5"	45.5"	47.5"	57.5"
Width	21.5"	21.5"	22"	22"	23.5"	23.5"	23.5"	24.5"
Height (ground to top of back of seat)	30.375" - 32.375"	34.75" - 36.75"	37" - 39"	37" - 39"	36.5" - 39.5"	36.5" - 39.5"	36.5" - 39.5"	42.5" - 45.5"
Height (ground to top of back of headrest)	N/A	N/A	N/A	N/A	N/A	N/A	N/A	50" - 56"
Ground to Top of Seat	19.25" - 21.25"	22" - 24"	23.5" - 25.5"	23.5" - 25.5"	23" - 26"	23" - 26"	23" - 26"	25" - 28"
Ground to Top of Deck	4.25"	4.125"	6"	6"	6.25"	6.25"	6.25"	6.5"
Deck to Top of Seat	15" - 17"	17.25" - 19.25"	18" - 20"	18" - 20"	17" - 20"	17" - 20"	17" - 20"	18.75" - 21.75"
Number of Seat Height Positions	3	4	3	3	4	4	4	4
Size of Increments	1"	5"	1"	1"	1"	1"	1"	1"
Front Axle to Rear Axle	28.625"	30"	30"	30"	30"	34.5"	35.5"	38.5"
Weight of Unit (assembled)	96 lbs	102 lbs	121 lbs	129 lbs	169 lbs	173 lbs	186 lbs	268 lbs
Front Half	29 lbs	33 lbs	44 lbs	52 lbs	46 lbs	48 lbs	63 lbs	68 lbs
Rear Half	26 lbs	25 lbs	28 lbs	28 lbs	48 lbs	48 lbs	48 lbs	68 lbs
Base	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A
Seat w/Arms	21 lbs	21 lbs	22 lbs	22 lbs	28 lbs	28 lbs	28 lbs	55 lbs
Batteries	19 lbs/pack	21 lbs/pack	25 lbs/pack	25 lbs/pack	22 lbs ea	22 lbs ea	22 lbs ea	37 lbs ea
Tires:	Solid	Solid	Solid, low profile	Solid, low profile	Solid, low profile	Solid, low profile	Solid, low profile	Pneumatic
Front	8"	8"	9"	9"	9"	9"	10"	13"
Rear	8"	8"	9"	9"	9"	N/A	N/A	13"
Casters	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A
Rear Anti-Tippers	2"	2"	2"	2"	2"	3"	3"	3"
Standard Seat:	Stadium Style	Stadium Style	Stadium Style	Stadium Style	Stadium Style	Stadium Style	Stadium Style	Captain's Seat
Slide Seat	N/A	N/A	N/A	N/A	Standard	Standard	Standard	N/A
Back Height - no headrest	16"	16"	16"	16"	16"	16"	16"	21" - 24"
Back Height with headrest	N/A	N/A	N/A	N/A	N/A	N/A	N/A	27.5" - 31.5"
Width x Depth	17" X 16"	17" X 16"	17" X 16"	17" X 16"	18" X 16"	18" X 16"	18" X 16"	20" X 17"
Color	Black Vinyl	Black Vinyl	Black Vinyl	Black Vinyl	Black Vinyl	Black Vinyl	Black Vinyl	Gray Vinyl
Optional Seats:	Stadium Style	Stadium Style	Stadium Style	Stadium Style	Stadium Style	Stadium Style	Stadium Style	N/A
Back Height - no headrest	16"	16"	16"	16"	16"	16"	16"	N/A
Back Height with headrest	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A
Width X Depth	18" X 16", 20" X 18"	18" X 16", 20" X 18"	18" X 16", 20" X 18"	18" X 16", 20" X 18"	20" X 18"	20" X 18"	20" X 18"	N/A
Color	Black Vinyl	Black Vinyl	Black	Black	Black	Black	Black	N/A
Optional Seats 2:	N/A	N/A	N/A	N/A	N/A	Captain's Highback	Captain's Highback	N/A
Back Height - no headrest	N/A	N/A	N/A	N/A	N/A	21"	21"	N/A
Back Height with headrest	N/A	N/A	N/A	N/A	N/A	29"	29"	N/A
Width	N/A	N/A	N/A	N/A	N/A	18"	18"	N/A
Depth	N/A	N/A	N/A	N/A	N/A	18"	18"	N/A
Color	N/A	N/A	N/A	N/A	N/A	Coffee/Sand, Charcoal/Gray	Coffee/Sand, Charcoal/Gray	N/A
Seat Options:								
Pneumatic Rise Seat Option	N/A	N/A	N/A	N/A	Yes	Yes	Yes	N/A
Electric Rise Seat Option	N/A	N/A	N/A	N/A	Yes	Yes	Yes	N/A
Footrest:	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A
Height adjustable	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A
Angle adjustable	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A
Warranty:								
Frame	Lifetime	Lifetime	Lifetime	Lifetime	Lifetime	Lifetime	Lifetime	Lifetime
Drive Train	2 Years	2 Years	2 Years	2 Years	2 Years	2 Years	2 Years	2 Years
Electronics ²	2 Years	2 Years	2 Years	2 Years	2 Years	2 Years	2 Years	2 Years

1. Battery range will vary due to rider weight, drive surface and terrain.
 2. Electronics warranty excludes batteries.
 3. The GP615 is sold separately from seat.
 4. Range tested at 200 lbs.

Literature is current at time of printing. Golden Technologies reserves the right to make changes to the product or literature at any time.

Power Chair Line Specifications

Roll



	ALANTÉ				COMPASS									
Model	GP204	GP215	GP605CC	GP605SS	GP600CC	GP600SS	GP615CC	GP615SS	GP620CC	GP620SS	GP600CC	GP600SS	GP615CC	GP615SS
Model Name	Alanté DX	Alanté LT	Compass Sport	Compass Sport	Compass	Compass	Compass TRO	Compass TRO	Compass HD	Compass HD	Compass	Compass	Compass TRO	Compass TRO
Medicare Code	K0823	K0816	K0823	K0822	K0823	K0822	K0836	K0835	K0825	K0824	K0823	K0822	K0836	K0835
Weight Capacity	300 lbs	300 lbs	300 lbs	300 lbs	300 lbs	300 lbs	300 lbs	300 lbs	450 lbs	450 lbs	300 lbs	300 lbs	300 lbs	300 lbs
Drive Wheels	Front	Front	Center	Center	Center	Center	Center	Center	Center	Center	Center	Center	Center	Center
Maximum Speed	3.5 mph	4.5 mph	4 mph	4 mph	5 mph	5 mph	4 mph	4 mph	4 mph	4 mph	4 mph	4 mph	4 mph	4 mph
Operating Range ¹	8 miles	6.5 miles	15 miles ⁴	15 miles ⁴	20 miles	20 miles	12 Miles	12 Miles	12 Miles	12 Miles	12 Miles	12 Miles	12 Miles	12 Miles
Ground Clearance	2.25"	2.25"	3.25"	3.25"	3.5"	3.5"	3.5"	3.5"	3.5"	3.5"	3.5"	3.5"	3.5"	3.5"
Turning Radius	24.5"	24.5"	19.5"	19.5"	20"	20"	20"	20"	20"	20"	20"	20"	20"	20"
Type Batteries	2-22AH	2-18AH	2-U1	2-U1	2-U1 or 2-NF22	2-U1 or 2-NF22	2-NF22	2-NF22	2-NF22	2-NF22	2-NF22	2-NF22	2-NF22	2-NF22
Colors	Red	Red	Red, Blue	Red, Blue	Red, Blue, Green, Onyx	Red, Blue, Green, Onyx	Red, Blue, Green, Onyx	Red, Blue, Green, Onyx	Red, Blue, Green, Onyx	Red, Blue, Green, Onyx	Red, Blue, Green, Onyx	Red, Blue, Green, Onyx	Red, Blue, Green, Onyx	Red, Blue, Green, Onyx
Freewheel Mode	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes
Electronic Speed Control	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes
Electro-Mechanical Brakes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes
Charger	Off Board	Off Board	Off Board	Off Board	Off Board	Off Board	Off Board	Off Board	Off Board	Off Board	Off Board	Off Board	Off Board	Off Board
Controller	PG VR-2 50 Amp	PG VR-2 50 Amp	Dynamic Shark 60 Amp Remote	Dynamic Shark 60 Amp Remote	Dynamic Shark 60 Amp Remote	Dynamic Shark 60 Amp Remote	Dynamic 75 Amp Remote	Dynamic 75 Amp Remote	Dynamic 75 Amp Remote	Dynamic 75 Amp Remote	Dynamic 75 Amp Remote	Dynamic 75 Amp Remote	Dynamic 75 Amp Remote	Dynamic 75 Amp Remote
Dimensions:														
Length	40"	40"	39"	39"	38.5"	38.5"	40"	40"	40"	40"	40"	40"	40"	40"
Width	24"	24"	24"	24"	24.5"	24.5"	25"	25"	25"	25"	25"	25"	25"	25"
Height (ground to top of back of seat)	34.5" - 36.5"	34.5" - 36.5"	40" - 42"	40" - 42"	38.25" - 42.75"	38.25" - 42.75"	38.25" - 42.75"	38.25" - 42.75"	38.25" - 42.75"	38.25" - 42.75"	38.25" - 42.75"	38.25" - 42.75"	38.25" - 42.75"	38.25" - 42.75"
Height (ground to top of back of headrest)	41.5" - 43.5"	41.5" - 43.5"	47" - 52.5"	47" - 52.5"	45.25" - 53.75"	45.25" - 53.75"	47.25" - 53.75"	47.25" - 53.75"	47.25" - 53.75"	47.25" - 53.75"	47.25" - 53.75"	47.25" - 53.75"	47.25" - 53.75"	47.25" - 53.75"
Ground to Top of Seat	20.5" - 22.5"	20.5" - 22.5"	21" - 23"	21" - 23"	19.75" - 24.25"	19.75" - 24.25"	20" - 23"	20" - 23"	21.25" - 24.25"	21.25" - 24.25"	20" - 23"	20" - 23"	21.25" - 24.25"	21.25" - 24.25"
Ground to Top of Deck	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A
Deck to Top of Seat	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A
Number of Seat Height Positions	3	3	3	3	10	10	5	5	5	5	5	5	5	5
Size of Increments	1"	1"	1"	1"	0.5"	0.5"	0.5"	0.5"	0.5"	0.5"	0.5"	0.5"	0.5"	0.5"
Front Axle to Rear Axle	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A
Weight of Unit (assembled)	133 lbs	130 lbs	182 lbs	179 lbs	189 lbs	186 lbs	229 lbs	222 lbs	229 lbs	223 lbs	186 lbs	186 lbs	229 lbs	222 lbs
Front Half	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A
Rear Half	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A
Base	67 lbs	67 lbs	92 lbs	92 lbs	92 lbs	92 lbs	101 lbs	101 lbs	101 lbs	101 lbs	92 lbs	92 lbs	101 lbs	101 lbs
Seat w/Arms	36 lbs	36 lbs	43 lbs	40 lbs	49 lbs	46 lbs	51 lbs	48 lbs	51 lbs	48 lbs	46 lbs	46 lbs	51 lbs	48 lbs
Batteries	30 lbs ea.	27 lbs pack	22 lbs ea.	22 lbs ea.	22 lbs ea.	22 lbs ea.	37 lbs ea.	37 lbs ea.	37 lbs ea.	37 lbs ea.	37 lbs ea.	37 lbs ea.	37 lbs ea.	37 lbs ea.
Tires:	Flat Free	Flat Free	Flat Free	Flat Free	Flat Free	Flat Free	Flat Free	Flat Free	Flat Free	Flat Free	Flat Free	Flat Free	Flat Free	Flat Free
Front	9.5"	9.5"	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A
Rear	N/A	N/A	10"	10"	10"	10"	10"	10"	10"	10"	10"	10"	10"	10"
Casters	6"	6"	6"	6"	6"	6"	6"	6"	6"	6"	6"	6"	6"	6"
Front Anti-Tippers	4"	4"	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A
Standard Seat:	Captain's Highback	Captain's Highback	Captain's Highback	Captain's Van/Pan	Captain's Highback	Captain's Van/Pan	Captain's Highback ¹	Tilt Seating ¹	Captain's Highback	Captain's Van/Pan	Captain's Highback	Captain's Van/Pan	Captain's Highback ¹	Tilt Seating ¹
Slide Seat	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A
Back Height - no headrest	18"	18"	21"	23"	21"	23"	21"	21"	21"	21"	21"	21"	21"	21"
Back Height with headrest	24.5" - 26.5"	24.5" - 26.5"	28.5" - 32.5"	30.5" - 32.5"	28.5" - 32.5"	30.5" - 32.5"	28.5" - 32.5"	N/A	28.5" - 32.5"	30.5" - 32.5"	28.5" - 32.5"	30.5" - 32.5"	28.	