

# L200

Deluxe high strength lightweight (Dual axle)

Improvement and more deluxe for the high strength lightweight wheelchair, the 18" seat depth provides more support and comfortable while seated. Detachable backrest provides more compact package for shipping to save transportation expenses. Various accessories choices make this chair versatile to meet the needs of active end users.















## Features

- High-strength lightweight wheelchair, the weight of wheelchair without front rigging is less than 34lb
- Dual axles for converting into hemi height wheelchair
- Easy for disassembly and assembly the backrests
- Back angle adjustable, 90°, 95°, 100°, and 105°
- Height adjustable flip up desk arm with padded armrests
- Swing away in-out footrests with heel loops
- 24" x 1-3/8" PU tired Mag rear wheel with quick release axles
- 8" x 1" castors with aluminum forks
- Nylon seat and backrest upholsteries
- Limited lifetime warranty on frame

## ■ Dimensions(cm/in)

(kg/lb)

												
	seat		arms		back		footrests		overall dimensions (w/o front rigging)			total weight (w/o front rigging)
adult	45/18	40/16	49/19 <sup>1/2</sup>	25-40/10-16	40-45/16-18	38/15	51/20	79/31	96/38	68/27	31/12 <sup>1/2</sup>	15.3/34
narrow	40/16	40/16	49/19 <sup>1/2</sup>	25-40/10-16	40-45/16-18	38/15	51/20	79/31	96/38	63/25	31/12 <sup>1/2</sup>	14.9/33