

N700

Full recliner

Designed for users who need extra shoulder and trunk support. The full recliner enables the user to comfortably sleep in the wheelchair eliminating the necessity of transferring to a bed. The back is reclinable up to 180°.



Features













- Full recliner mechanism operated by and maglock, infinite angle adjusted easily
- Detachable desk arms (Full arms optional)
- Elevating leg rests
- Deluxe Nylon seat and long back upholsteries
- Headrest on extended tubes (pillow optional)
- 24" x 1" solid tired Mag wheel with chromed steel handrims
- 8" x 1" solid castors
- Anti-tippers
- Limited lifetime warranty on frame



Optional full arms

Dimensions(cm/in)

(kg/lb)

												
	seat			arms	back	footrests		overall dimensions (w/o front rigging)			total weight (w/o front rigging)	
wide	51/20	43/17	48/19	25/10	86/34	40/16	53/21	119/47	135/53	71/28	30/12	23/50.5
adult	45/18	43/17	48/19	25/10	86/34	40/16	53/21	119/47	135/53	66/26	30/12	22.2/49
narrow	40/16	43/17	48/19	25/10	86/34	40/16	53/21	119/47	135/53	61/24	30/12	21.5/47.5