

The **Invacare® Solara® 3G tilt-in-space wheelchair** is now available in specific configurations with some of the most popular options to take the guesswork out of configuring a complex wheelchair. With a standard weight capacity of 300 lb., back angle adjustment and seat width and depth growth, you can easily customize the **Solara 3G Wheelchair's** fit to match a specific individual's needs.

INVACARE® SOLARA® 3G WHEELCHAIR SmartSelect™ PROGRAM

Adjustability. Reliability. Serviceability.


Invacare® Solara® 3G Wheelchair SmartSelect Features:

- 1 **Color coded, rust-resistant zinc-plated hardware** make adjustments easier. The silver hardware indicates a point where an adjustment can occur while black bolts are static.
- 2 Each chair has **4" of simple seat width** and ample depth adjustability to accommodate a change in patient need.
- 3 The **seat pan system** is easily converted from flush mount to drop mount or vice versa utilizing the same hardware.
- 4 The unique **tilt assist feature** on all Solara 3G wheelchair makes tilting easier by reducing the effort required to tilt an occupied Solara 3G wheelchair. With three settings, front, middle and rear, the simple Tilt Assist adjustments help to correct for uneven weight distribution.
- 5 A **fully welded steel frame** provides a reliable tilt-in-space wheelchair for your facility.
- 6 **Transport Brackets (TRBKTS)** have been designed to allow for more access to the rear tie down points without additional steps needed.


For more information about chair specifications, contact your Invacare Sales Representative or visit, www.invacare.com/Rehab



Color coded hardware



Width & depth adjustability



Seat pan system



Transport brackets (TRBKTS)



Tilt assist feature

All frame color in photography are used to show chair features. Limited color selection available through Solara 3G Wheelchair SmartSelect programs.

©2012 Invacare Corporation. All rights reserved. Trademarks are identified by the symbols ™ and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Specifications are subject to change without notification. Form No. 12-286



Yes, you can.®